

effimed<sup>+</sup>



BIOMEDICAL FREEZERS

-25°C

# BIOMEDICAL FREEZERS -25°C

Efficold, with more than 30 years of experience as a manufacturer of cooling equipment and innovative solutions for the biomedical sector, already has divisions targeting the food, beverage, hospitality, supermarket, ice cream and special projects sectors. We also manufacture surgical masks and FFP2 masks.

We are part of the Onnera Group, a leading medical equipment manufacturing group. ONNERA GROUP is an international business group dedicated to providing equipment solutions for catering laundry, refrigeration applications and medical manufacturing.

FAGOR, EFFICOLD, ONNERA CONTRACT, EDENOX, DOMUS, ASBER, DANUBE and PRIMER are leaders in different sectors, markets and countries. A global group that markets its products on five continents and now also includes EFFIMED.



## LFV-100



Model	
-25° Biomedical Freezer	19086353
Model	LFV-100
Cabinet Type	VERTICAL
Capacity (Lt)	89
Internal Dimensions (W x D x H) (mm)	500 x 485 x 532
External Dimensions (W x D x H) (mm)	677 x 649 x 948
Packaging Dimensions (W x D x H) (mm)	765 x 735 x 935
Net Weight (Kg)	54

Performance	
Temperature Range	-10°C » -25°C
Ambient Temperature	25°C / 60% HR
Climatic Class	N
Controler	Microprocesator
Monitor	LCD Display

Refrigeration	
Compressors	1pz
Cooling Method	Direct
Defrost Mode	Manual
Refrigerant	R600a
Insulation Material/Thickness	100

Construction	
External Material	Powder coated material

Internal Material and Shelves	Aluminum plate with spraying
Shelves	2/3 Drawers
Outer Door Lock	Yes
Access Port	1
Weels and Caster	4/2
DATA/Interval/Recording Time	NO

Alarms	
Temperature Range	Yes
Electrical	Yes
System	Yes

Electrical	
Power Supply/Frequency	220V / 50Hz
Energy Consupion Kw/h 24 hrs	0.75