

effimed<sup>+</sup>



BIOMEDICAL FREEZERS

-25°C

# BIOMEDICAL FREEZERS -25°C

Efficold, with more than 30 years of experience as a manufacturer of cooling equipment and innovative solutions for the biomedical sector, already has divisions targeting the food, beverage, hospitality, supermarket, ice cream and special projects sectors. We also manufacture surgical masks and FFP2 masks.

We are part of the Onnera Group, a leading medical equipment manufacturing group. ONNERA GROUP is an international business group dedicated to providing equipment solutions for catering laundry, refrigeration applications and medical manufacturing.

FAGOR, EFFICOLD, ONNERA CONTRACT, EDENOX, DOMUS, ASBER, DANUBE and PRIMER are leaders in different sectors, markets and countries. A global group that markets its products on five continents and now also includes EFFIMED.



## LFV-450





## Model

-25? Biomedical Freezer	19084875
Model	LFV-450
Cabinet Type	VERTICAL
Capacity (Lt)	449
Internal Dimensions (W x D x H) (mm)	(650 x 570 x 627 ) x 2
External Dimensions (W x D x H) (mm)	810 x 735 x 1955
Packaging Dimensions (W x D x H) (mm)	895 x 820 x 2045
Net Weight (Kg)	133

## Performance

Temperature Range	-10°C » -25°C
Ambient Temperature	25°C / 60% HR
Climatic Class	N
Controler	Microprocesator
Monitor	LCD Display

## Refrigeration

Compressors	2pz
Cooling Method	Direct
Defrost Mode	Manual
Refrigerant	R600a
Insulation Material/Thickness	80

## Construction

External Material	Powder coated material
-------------------	------------------------

Internal Material and Shelves	Aluminum plate with spraying
Shelves	5/12 Drawers
Outer Door Lock	Yes
Access Port	2
Weels and Caster	4/2
DATA/Interval/Recording Time	USB/Record every 10 minute / 2 years

Alarms	
Temperature Range	Yes
Electrical	Yes
System	Yes

Electrical	
Power Supply/Frequency	220V / 50Hz
Energy Consupcion Kw/h 24 hrs	2.24