

effimed 



ULTRA FREEZERS


-40°C

ULTRA FREEZERS -40°C

Efficold, with more than 30 years of experience as a manufacturer of cooling equipment and innovative solutions for the biomedical sector, already has divisions targeting the food, beverage, hospitality, supermarket, ice cream and special projects sectors. We also manufacture surgical masks and FFP2 masks.

We are part of the Onnera Group, a leading medical equipment manufacturing group. ONNERA GROUP is an international business group dedicated to providing equipment solutions for catering laundry, refrigeration applications and medical manufacturing.

FAGOR, EFFICOLD, ONNERA CONTRACT, EDENOX, DOMUS, ASBER, DANUBE and PRIMER are leaders in different sectors, markets and countries. A global group that markets its products on five continents and now also includes EFFIMED.



MFC-300





Model

Ultra Freezer -40°C	19085274
Model	MFC-300
Cabinet Type	HORIZONTAL
Volume (Lt)	268
Capacity (Aluminium Racks 2" Boxes)	210
Internal Dimensions (W x D x H)	1019 x 465 x 651
External Dimensions (W x D x H)	1245 x 775 x 929
Packaging Dimensions (mm)	1349 x 815 x 975
Net Weight (Kg)	77

Performance

Temperature Setting Range	-10°C » -40°C
External Temperature	25°C / 60% HR
Climatic Class	N
Controler	Microprocesator
Display type	LED
Compressors	1pz
Cooling Method	Direct
Defrost Mode	Manual
Refrigerant	R290
Insulation Material/Thickness	110mm

Construction

External Material	SPCC epoxy coating
-------------------	--------------------

Internal Material and Shelves	Stainless steel
Shelves/Drawers	1
Number of Inner Doors	0
Outer Door Lock	Yes
Access Port	1
Weels and Caster	4/2
DATA/Interval/	USB/10 minutes /2 years
Battery back-up	Yes

Alarms

System Alarm	Yes
Power supply alarm	Yes
System Alarm	Yes

Electrical

Power Supply/Frequency	220V / 50Hz
Power Consumption(KWh/24h)	3.1
110V-220V 60 Hz Version	Optional

Sound level

Noise Level Db	43
----------------	----