

effimed⁺



ULTRA FREEZERS

-40°C

ULTRA FREEZERS -40°C

Efficold, with more than 30 years of experience as a manufacturer of cooling equipment and innovative solutions for the biomedical sector, already has divisions targeting the food, beverage, hospitality, supermarket, ice cream and special projects sectors. We also manufacture surgical masks and FFP2 masks.

We are part of the Onnera Group, a leading medical equipment manufacturing group. ONNERA GROUP is an international business group dedicated to providing equipment solutions for catering laundry, refrigeration applications and medical manufacturing.

FAGOR, EFFICOLD, ONNERA CONTRACT, EDENOX, DOMUS, ASBER, DANUBE and PRIMER are leaders in different sectors, markets and countries. A global group that markets its products on five continents and now also includes EFFIMED.



MFC-400



Model		
Ultra Freezer -40°C		19084876
Model		MFC-400
Cabinet Type		HORIZONTAL
Volume (Lt)		358
Capacity (Aluminium Racks 2" Boxes)		270
Internal Dimensions (W x D x H)		1308 x 465 x 651
External Dimensions (W x D x H)		1534 x 775 x 929
Packaging Dimensions (mm)		1783 x 815 x 975
Net Weight (Kg)		98

Performance		
Temperature Setting Range		-10°C » -40°C
External Temperature		25°C / 60% HR
Climatic Class		N
Controler		Microprocesator
Display type		LED
Compressors		1pz
Cooling Method		Direct
Defrost Mode		Manual
Refrigerant		R290
Insulation Material/Thickness		110mm

Construction		
External Material		SPCC epoxy coating

Internal Material and Shelves	Stainless steel
Shelves/Drawers	1
Number of Inner Doors	0
Outer Door Lock	Yes
Access Port	1
Weels and Caster	4/2
DATA/Interval/	USB / 10 minutes / 2 years
RS 485 Port	No

Alarms

System Alarm	Yes
Power supply alarm	Yes
System Alarm	Yes

Electrical

Power Supply/Frequency	220V / 50Hz
Power Consumption(KWh/24h)	3.88
110V-220V 60 Hz Version	Optional

Sound level

Noise Level Db	43
----------------	----