

effimed⁺



ULTRA FREEZERS

-40°C

ULTRA FREEZERS -40°C

Efficold, with more than 30 years of experience as a manufacturer of cooling equipment and innovative solutions for the biomedical sector, already has divisions targeting the food, beverage, hospitality, supermarket, ice cream and special projects sectors. We also manufacture surgical masks and FFP2 masks.

We are part of the Onnera Group, a leading medical equipment manufacturing group. ONNERA GROUP is an international business group dedicated to providing equipment solutions for catering laundry, refrigeration applications and medical manufacturing.

FAGOR, EFFICOLD, ONNERA CONTRACT, EDENOX, DOMUS, ASBER, DANUBE and PRIMER are leaders in different sectors, markets and countries. A global group that markets its products on five continents and now also includes EFFIMED.



MFV-300



Model	
Ultra Freezer -40°C	19086356
Model	MFV-300
Cabinet Type	VERTICAL
Volume (Lt)	268
Capacity (Aluminium Racks 2" Boxes)	-
Internal Dimensions (W x D x H)	500 x 460 x 1235
External Dimensions (W x D x H)	700 x 640 x 1792
Packaging Dimensions (mm)	757 x 717x 1885
Net Weight (Kg)	94

Performance	
Temperature Setting Range	-10°C » -40°C
External Temperature	25°C / 60% HR
Climatic Class	N
Controler	Microprocesator
Display type	LCD screen
Compressors	1pz
Cooling Method	Direct
Defrost Mode	Manual
Refrigerant	R290
Insulation Material/Thickness	100mm

Construction	
External Material	Galvanized

Internal Material and Shelves	aluminium
Shelves/Drawers	7/7
Number of Inner Doors	0
Outer Door Lock	Yes
Access Port	1
Weels and Caster	4/2
DATA/Interval/	Screen data
Battery back-up	No

Alarms

System Alarm	No
Power supply alarm	Yes
System Alarm	Yes

Electrical

Power Supply/Frequency	220V / 50Hz
Power Consumption(KWh/24h)	1.69
110V-220V 60 Hz Version	Optional

Sound level

Noise Level Db	40
----------------	----