

effimed⁺



ULTRA FREEZERS

-40°C

ULTRA FREEZERS -40°C

Efficold, with more than 30 years of experience as a manufacturer of cooling equipment and innovative solutions for the biomedical sector, already has divisions targeting the food, beverage, hospitality, supermarket, ice cream and special projects sectors. We also manufacture surgical masks and FFP2 masks.

We are part of the Onnera Group, a leading medical equipment manufacturing group. ONNERA GROUP is an international business group dedicated to providing equipment solutions for catering laundry, refrigeration applications and medical manufacturing.

FAGOR, EFFICOLD, ONNERA CONTRACT, EDENOX, DOMUS, ASBER, DANUBE and PRIMER are leaders in different sectors, markets and countries. A global group that markets its products on five continents and now also includes EFFIMED.



MFV-450



Model

Model	MFV-450
Cabinet Type	VERTICAL
Volume (Lt)	435
Capacity (Aluminium Racks 2" Boxes)	-
Internal Dimensions (W x D x H)	650 x 547 x 1305
External Dimensions (W x D x H)	850 x 740 x 1935
Packaging Dimensions (mm)	990 x 855 x 2095
Net Weight (Kg)	175

Performance

Temperature Setting Range	-10°C » -40°C
External Temperature	25°C / 60% HR
Climatic Class	N
Controler	Microprocesator
Display type	Digital Display
Compressors	1pz
Cooling Method	Direct
Defrost Mode	Manual
Refrigerant	R507
Insulation Material/Thickness	100mm

Construction

External Material	Carbon steel spraying
Internal Material and Shelves	Stainless steel

Shelves/Drawers	7/14
Number of Inner Doors	0
Outer Door Lock	Yes
Access Port	1
Weels and Caster	4/2
DATA/Interval/	USB / 10 minutes / 2 years
Battery back-up	Yes

Alarms

System Alarm	Yes
Power supply alarm	Yes
System Alarm	Yes

Electrical

Power Supply/Frequency	220V / 50Hz
Power Consumption(KWh/24h)	4.68
110V-220V 60 Hz Version	Optional

Sound level

Noise Level Db	44
----------------	----