

effimed⁺



ULTRA FREEZERS

-40°C

ULTRA FREEZERS -40°C

Efficold, with more than 30 years of experience as a manufacturer of cooling equipment and innovative solutions for the biomedical sector, already has divisions targeting the food, beverage, hospitality, supermarket, ice cream and special projects sectors. We also manufacture surgical masks and FFP2 masks.

We are part of the Onnera Group, a leading medical equipment manufacturing group. ONNERA GROUP is an international business group dedicated to providing equipment solutions for catering laundry, refrigeration applications and medical manufacturing.

FAGOR, EFFICOLD, ONNERA CONTRACT, EDENOX, DOMUS, ASBER, DANUBE and PRIMER are leaders in different sectors, markets and countries. A global group that markets its products on five continents and now also includes EFFIMED.



MFV-800



Model	
Ultra Freezer -40°C	19084879
Model	MFV-800
Cabinet Type	VERTICAL
Volume (Lt)	775
Capacity (Aluminium Racks 2" Boxes)	600
Internal Dimensions (W x D x H)	865 x 696 x 1286
External Dimensions (W x D x H)	1205 x 1025 x 1955
Packaging Dimensions (mm)	1320 x 1155 x 2171
Net Weight (Kg)	258

Performance	
Temperature Setting Range	-10°C » -40°C
External Temperature	25°C / 60% HR
Climatic Class	N
Controler	Microprocesator
Display type	Digital Display
Compressors	2pz
Cooling Method	Direct
Defrost Mode	Manual
Refrigerant	Mixture
Insulation Material/Thickness	130mm

Construction	
External Material	Steel plate with spraying

Internal Material and Shelves	AISI 304
Shelves/Drawers	3/0
Number of Inner Doors	4
Outer Door Lock	Yes
Access Port	3
Weels and Caster	4/2
DATA/Interval/	USB / 10 minutes / 2 years
Battery back-up	Yes

Alarms

System Alarm	Yes
Power supply alarm	Yes
System Alarm	Yes

Electrical

Power Supply/Frequency	220V / 50Hz
Power Consumption(KWh/24h)	7.02
110V-220V 60 Hz Version	Optional

Sound level

Noise Level Db	45
----------------	----