

effimed 



ULTRA FREEZERS

-440°C

# ULTRA FREEZERS -40°C

Efficold, with more than 30 years of experience as a manufacturer of cooling equipment and innovative solutions for the biomedical sector, already has divisions targeting the food, beverage, hospitality, supermarket, ice cream and special projects sectors. We also manufacture surgical masks and FFP2 masks.

We are part of the Onnera Group, a leading medical equipment manufacturing group. ONNERA GROUP is an international business group dedicated to providing equipment solutions for catering laundry, refrigeration applications and medical manufacturing.

FAGOR, EFFICOLD, ONNERA CONTRACT, EDENOX, DOMUS, ASBER, DANUBE and PRIMER are leaders in different sectors, markets and countries. A global group that markets its products on five continents and now also includes EFFIMED.



## MFV-90





## Model

Ultra Freezer -40°C	19086355
Model	MFV-90
Cabinet Type	VERTICAL
Volume (Lt)	90
Capacity (Aluminium Racks 2" Boxes)	-
Internal Dimensions (W x D x H)	500 x 485 x 532
External Dimensions (W x D x H)	677 x 649 x 948
Packaging Dimensions (mm)	765 x 735 x 935
Net Weight (Kg)	54

## Performance

Temperature Setting Range	-10°C » -40°C
External Temperature	25°C / 60% HR
Climatic Class	N
Controler	Microprocesator
Display type	Digital Display
Compressors	1pz
Cooling Method	Direct
Defrost Mode	Manual
Refrigerant	R290
Insulation Material/Thickness	100mm

## Construction

External Material	Powder coater material
-------------------	------------------------

Internal Material and Shelves	Aluminium plate with spray
Shelves/Drawers	2/3
Number of Inner Doors	0
Outer Door Lock	Yes
Access Port	1
Weels and Caster	4/2
DATA/Interval/	Screen data
Battery back-up	No

### Alarms

System Alarm	Yes
Power supply alarm	No
System Alarm	No

### Electrical

Power Supply/Frequency	220V / 50Hz
Power Consumption(KWh/24h)	1.69
110V-220V 60 Hz Version	Optional

### Sound level

Noise Level Db	40
----------------	----