

effimed⁺



COMBINED FREEZING / REFRIGERATION

COMBIS

ULTRA FREEZING-REFRIGERATION

-10°C ~ -40°C // 2°C ~ 8°C

Efficold, with more than 30 years of experience as a manufacturer of cooling equipment and innovative solutions for the biomedical sector, already has divisions targeting the food, beverage, hospitality, supermarket, ice cream and special projects sectors. We also manufacture surgical masks and FFP2 masks.

We are part of the Onnera Group, a leading medical equipment manufacturing group. ONNERA GROUP is an international business group dedicated to providing equipment solutions for catering laundry, refrigeration applications and medical manufacturing.

FAGOR, EFFICOLD, ONNERA CONTRACT, EDENOX, DOMUS, ASBER, DANUBE and PRIMER are leaders in different sectors, markets and countries. A global group that markets its products on five continents and now also includes EFFIMED.



PC-450



Model	
Combi Medical Cooler Ultra Freezer -40°C	19084872
Model	PC-450
Cabinet Type	VERTICAL
Volume (C/F) Litres	R: 223 - F: 223
Internal Dimensions (W x D x H) mm	Refrigerator: 650 x 570 x 627
External Dimensions (W x D x H) mm	810 x 745 x 1955
Packaging Dimensions (W x D x H) mm	895 x 820 x 2035
Net Weight (Kg)	146

Performance	
Temperature Setting Range	4°C » 8°C / 10°C » -40°C
External Temperature	25°C / 60%HR
Climatic Class	N
Controler	Microprocesator
Display type	Digital Display

Refrigeration	
Compresor	2pz
Cooling Method	R: Forced air cooling - F: Direct cooling
Defrost Mode	R: Automatic - F: Manual
Refrigerant	R: R600a - F:290
Insulation Thickness (mm)	R: 80 - F: 80

Construction	
External Material	Powder coated material

Internal Material and Shelves	Aluminum plate with spraying
Shelves	R: 3 (coated steel wired shelf) - F: 6 (ABS)
Outer Door Lock	Yes
Heated Sealed Door	Yes
Access Port	2
Weels and Caster	4/2
DATA/Interval/Recording Time	USB / Record every 10 minute / 2 years
Battery	Yes

Alarms

Temperature Setting Range	Yes
Electrical	Yes
System	Yes

Electrical

Voltage / frequency	220V / 50Hz
24h energy consumption	4.95

Sound level

Noise Level (Db)	43
------------------	----