



effimed 



ULTRA LOW TEMPERATURE FREEZERS

 **-86°C** 

ULTRA LOW TEMPERATURE FREEZERS -86°C

Efficold, with more than 30 years of experience as a manufacturer of cooling equipment and innovative solutions for the biomedical sector, already has divisions targeting the food, beverage, hospitality, supermarket, ice cream and special projects sectors. We also manufacture surgical masks and FFP2 masks.

We are part of the Onnera Group, a leading medical equipment manufacturing group. ONNERA GROUP is an international business group dedicated to providing equipment solutions for catering laundry, refrigeration applications and medical manufacturing.

FAGOR, EFFICOLD, ONNERA CONTRACT, EDENOX, DOMUS, ASBER, DANUBE and PRIMER are leaders in different sectors, markets and countries. A global group that markets its products on five continents and now also includes EFFIMED.



UFC-350





Model

Ultra Freezer -86°C	19084880
Model	UFC-350
Cabinet Type	HORIZONTAL
Volume (Lt)	326
Capacity (Aluminium Racks 2" Boxes)	216
Internal Dimensions (W x D x H)	1200 x 470 x 570
External Dimensions (W x D x H)	2005 x 890 x 1009
Packaging Dimensions (mm)	2193 x 1003 x 1256
Net Weight	280

Performance

Temperature Setting Range	-40°C » -86°C
Ambient Temperature	25°C / 60% HR
Climatic Class	N
Controler	Microprocesator
Display type	Digital Display

Cooling

Compressors	1pz
Cooling Method	Direct
Defrost Mode	Manual
Refrigerant	Mixture

Construction

Insulation Material/Thickness	130 mm
-------------------------------	--------

External Material	Steel plate with spraying
Internal Material and Shelves	AISI 304
Shelves	0
Number of Inner Doors	0
Outer Door Lock	Yes
Vacuum Release Port	No
Heated Seal Door	No
Access Port	1
Weels and Caster	Yes
DATA/Interval/Recording Time	Printer
Batery/Hrs	Yes

Alarms

Stock alarm	-
-------------	---

Electrical

Power Supply / Frequency	220V / 50Hz
--------------------------	-------------

Sound level

Noise Level Db	51
----------------	----