

effimed<sup>+</sup>



ULTRA LOW TEMPERATURE FREEZERS

-86°C

# ULTRA LOW TEMPERATURE FREEZERS -86°C

Efficold, with more than 30 years of experience as a manufacturer of cooling equipment and innovative solutions for the biomedical sector, already has divisions targeting the food, beverage, hospitality, supermarket, ice cream and special projects sectors. We also manufacture surgical masks and FFP2 masks.

We are part of the Onnera Group, a leading medical equipment manufacturing group. ONNERA GROUP is an international business group dedicated to providing equipment solutions for catering laundry, refrigeration applications and medical manufacturing.

FAGOR, EFFICOLD, ONNERA CONTRACT, EDENOX, DOMUS, ASBER, DANUBE and PRIMER are leaders in different sectors, markets and countries. A global group that markets its products on five continents and now also includes EFFIMED.



## UFC-350



Model		
Ultra Freezer -86°C		19084880
Model		UFC-350
Cabinet Type		HORIZONTAL
Volume (Lt)		326
Capacity (Aluminium Racks 2" Boxes)		216
Internal Dimensions (W x D x H)		1200 x 470 x 570
External Dimensions (W x D x H)		2005 x 890 x 1009
Packaging Dimensions (mm)		2193 x 1003 x 1256
Net Weight		280

Performance		
Temperature Setting Range		-40°C » -86°C
Ambient Temperature		25°C / 60% HR
Climatic Class		N
Controler		Microprocesator
Display type		Digital Display

Cooling		
Compressors		1pz
Cooling Method		Direct
Defrost Mode		Manual
Refrigerant		Mixture

Construction		
Insulation Material/Thickness		130 mm

External Material	Steel plate with spraying
Internal Material and Shelves	AISI 304
Shelves	0
Number of Inner Doors	0
Outer Door Lock	Yes
Vacuum Release Port	No
Heated Seal Door	No
Access Port	1
Weels and Caster	Yes
DATA/Interval/Recording Time	Printer
Batery/Hrs	Yes

## Alarms

Stock alarm	-
-------------	---

## Electrical

Power Supply / Frequency	220V / 50Hz
--------------------------	-------------

## Sound level

Noise Level Db	51
----------------	----