

effimed<sup>+</sup>



ULTRA LOW TEMPERATURE FREEZERS

-86°C

# ULTRA LOW TEMPERATURE FREEZERS -86°C

Efficold, with more than 30 years of experience as a manufacturer of cooling equipment and innovative solutions for the biomedical sector, already has divisions targeting the food, beverage, hospitality, supermarket, ice cream and special projects sectors. We also manufacture surgical masks and FFP2 masks.

We are part of the Onnera Group, a leading medical equipment manufacturing group. ONNERA GROUP is an international business group dedicated to providing equipment solutions for catering laundry, refrigeration applications and medical manufacturing.

FAGOR, EFFICOLD, ONNERA CONTRACT, EDENOX, DOMUS, ASBER, DANUBE and PRIMER are leaders in different sectors, markets and countries. A global group that markets its products on five continents and now also includes EFFIMED.



## UFV-400



Model		
Ultra Freezer -86°C		19084882
Model		UFV-400
Cabinet Type		VERTICAL
Volume (Lt)		400
Capacity (Aluminium Racks 2" Boxes)		300
Internal Dimensions (W x D x H)		440 x 696 x 1266
External Dimensions (W x D x H)		785 x 1041 x 1947
Packaging Dimensions (mm)		890 x 1165 x 2145
Net Weight		237

Performance		
Temperature Setting Range		-40°C » -86°C
Ambient Temperature		25°C / 60% HR
Climatic Class		N
Controler		Microprocesator
Display type		Touch screen

Cooling		
Compressors		1pz
Cooling Method		Direct
Defrost Mode		Manual
Refrigerant		Combined

Construction		
Insulation Material/Thickness		130 mm

External Material	Prepainted steel
Internal Material and Shelves	AISI 304
Shelves	3
Number of Inner Doors	2
Outer Door Lock	Yes
Vacuum Release Port	Yes
Heated Seal Door	Yes
Access Port	2
Weels and Caster	Yes
DATA/Interval/Recording Time	USB -
Batery/Hrs	Yes

## Alarms

Stock alarm	Yes
-------------	-----

## Electrical

Power Supply / Frequency	220V / 50Hz
--------------------------	-------------

## Sound level

Noise Level Db	50
----------------	----