

effimed 



ULTRA LOW TEMPERATURE FREEZERS

-86°C

ULTRA LOW TEMPERATURE FREEZERS -86°C

Efficold, with more than 30 years of experience as a manufacturer of cooling equipment and innovative solutions for the biomedical sector, already has divisions targeting the food, beverage, hospitality, supermarket, ice cream and special projects sectors. We also manufacture surgical masks and FFP2 masks.

We are part of the Onnera Group, a leading medical equipment manufacturing group. ONNERA GROUP is an international business group dedicated to providing equipment solutions for catering laundry, refrigeration applications and medical manufacturing.

FAGOR, EFFICOLD, ONNERA CONTRACT, EDENOX, DOMUS, ASBER, DANUBE and PRIMER are leaders in different sectors, markets and countries. A global group that markets its products on five continents and now also includes EFFIMED.



UFV-800





Model

Ultra Freezer -86°C	19084884
Model	UFV-800
Cabinet Type	VERTICAL
Volume (Lt)	775
Capacity (Aluminium Racks 2" Boxes)	600
Internal Dimensions (W x D x H)	865 x 696 x 1286
External Dimensions (W x D x H)	1205 x 1025 x 1955
Packaging Dimensions (mm)	1320 x 1155 x 2171
Net Weight	365

Performance

Temperature Setting Range	-40°C » -86°C
Ambient Temperature	25°C / 60% HR
Climatic Class	N
Controler	Microprocesator
Display type	Touch screen

Cooling

Compressors	2pz
Cooling Method	Direct
Defrost Mode	Manual
Refrigerant	Combined

Construction

Insulation Material/Thickness	130 mm
-------------------------------	--------

External Material	Prepainted steel
Internal Material and Shelves	AISI 304
Shelves	3
Number of Inner Doors	4
Outer Door Lock	Yes
Vacuum Release Port	Yes
Heated Seal Door	Yes
Access Port	3
Weels and Caster	Yes
DATA/Interval/Recording Time	USB -
Batery/Hrs	Yes

Alarms

Stock alarm	Yes
-------------	-----

Electrical

Power Supply / Frequency	220V / 50Hz
--------------------------	-------------

Sound level

Noise Level Db	50
----------------	----