

effimed⁺



ULTRA LOW TEMPERATURE FREEZERS

-86°C

ULTRA LOW TEMPERATURE FREEZERS -86°C

Efficold, with more than 30 years of experience as a manufacturer of cooling equipment and innovative solutions for the biomedical sector, already has divisions targeting the food, beverage, hospitality, supermarket, ice cream and special projects sectors. We also manufacture surgical masks and FFP2 masks.

We are part of the Onnera Group, a leading medical equipment manufacturing group. ONNERA GROUP is an international business group dedicated to providing equipment solutions for catering laundry, refrigeration applications and medical manufacturing.

FAGOR, EFFICOLD, ONNERA CONTRACT, EDENOX, DOMUS, ASBER, DANUBE and PRIMER are leaders in different sectors, markets and countries. A global group that markets its products on five continents and now also includes EFFIMED.



UFV-PORTABLE



Model

Ultra Freezer -86°C	19085162
Model	UFV-PORTABLE
Cabinet Type	PORTABLE
Volume (Lt)	1.8
Capacity (Aluminium Racks 2" Boxes)	0
Internal Dimensions (W x D x H)	152 x 133 x 87
External Dimensions (W x D x H)	245 x 282 x 496
Packaging Dimensions (mm)	305 x 295 x 570
Net Weight	12

Performance

Temperature Setting Range	-40°C » -86°C
Ambient Temperature	25°C / 60%HR
Climatic Class	N
Controler	Microprocesator
Display type	LCD screen

Cooling

Compressors	Peltier
Cooling Method	Direct
Defrost Mode	Manual
Refrigerant	-

Construction

Insulation Material/Thickness	-
-------------------------------	---

External Material	Prepainted steel
Internal Material and Shelves	EVA
Shelves	0
Number of Inner Doors	0
Outer Door Lock	Yes
Vacuum Release Port	No
Heated Seal Door	No
Access Port	0
Weels and Caster	Yes
DATA/Interval/Recording Time	Data in screen
Batery/Hrs	No

Alarms

Stock alarm	-
-------------	---

Electrical

Power Supply / Frequency	220V / 50Hz
--------------------------	-------------

Sound level

Noise Level Db	45
----------------	----